

Get confidence from your customer with successful litter



Complete extender range for dog semen

CANIXcell

helps maintain sperm cells viable for 24 h



Easy to use

1 ejaculate = 1 bottle

Canine semen cooling extender:

CANIXCell: ref 021729



CANIFreeze

Enhanced formula for frozen dog semen



Canine semen freezing extender:

CANIFreeze: ref 021728




Extenders ready to use ⇨ nothing to add

Improved formulae ⇨ contains egg yolk plasma

Sterile

Proven results from the veterinary university of Nantes (F) (ESDAR 2008)

WORLD LEADER IN REPRODUCTION BIOTECHNOLOGIES



www.imv-technologies.com